

## CALIBRATION RECORD OF DISINTEGRATION TEST APPARATUS

Page No.: 1 of 1

<b>Name</b>	<b>Disintegration Test Apparatus</b>	<b>Location</b>	<b>IPQA Lab.</b>
<b>Make / Model</b>		<b>Frequency</b>	<b>Monthly</b>
<b>ID No.</b>		<b>Calibration done on</b>	
<b>Serial No.</b>		<b>Next Due Date</b>	

### 1) CALIBRATION OF THERMOSTAT:

Calibrated Thermometer ID No.: \_\_\_\_\_

Sr. No.	Time (Min)	Observation				Acceptance Limits (37°C ± 2°C)
		LHS Beaker (°C)		RHS Beaker (°C)		
		Thermometer	Display	Thermometer	Display	
1	05					35°C to 39°C
2	15					35°C to 39°C
3	30					35°C to 39°C
4	60					35°C to 39°C
Remarks:				Done By (Sign./Date):		

### 2) CALIBRATION OF UP AND DOWN CYCLES PER MINUTE:

Calibrated Stop watch ID No.: \_\_\_\_\_

Sr. No.	Observation (Cycles per minute)	Limits (Cycles per minute)	Done By (Sign./Date)
1		29 to 32	
2		29 to 32	
3		29 to 32	
Remarks (If Any):			

**Remarks: The Disintegration Test Apparatus is calibrated and is Qualified / Not Qualified for use.**

<b>Done By</b>	<b>Checked By</b>	<b>Verified By (QA)</b>
<b>(Sign./Date)</b>	<b>(Sign./Date)</b>	<b>(Sign./Date)</b>