

VISUAL OBSERVATION OF CONTROL SAMPLE FOR CAPSULE

Page 1 of 1

Product Brand/Generic Name: _____ Batch No.: _____

MFG Date: _____ Expiry Date: _____

Observation Year	Defects						Sign and date	Destroyed by	Remarks
	Appearance (free from mashed, bent, double cap, scrape, improper print (Ok/Not OK)	Dented / color variation / discoloration / swelling/split capsule (Observed / Not observed)	Locking defect / loose cap / Body or cap broken /crack capsule / black spot on body or cap (Observed/ Not observed)	Telescopic / foreign particle/ peeling / punched ends (Observed/ Not observed)	Printing on Carton / Blister / Strip (Ok/Not OK)	Others			
1 st Year _____									
2 nd Year _____									
3 rd Year _____									
4 th Year _____									
Other									